

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High		Low	
Mean	27.18	Mean	28.04	Mean	25.76
SD	0.37	SD	0.4	SD	0.26
CV	1.35	CV	1.43	CV	0.99
Day-Run	Value	Day-Run	Value	Day-Run	Value
2 - 1	27.6	2 - 1	28.6	2 - 1	26
3 - 1	27	3 - 1	28.5	3 - 1	26
4 - 1	27.5	4 - 1	28.3	4 - 1	25.9
5 - 1	26.9	5 - 1	28.1	5 - 1	25.7
6 - 1	27.6	6 - 1	28.2	6 - 1	25.5
9 - 1	27.5	9 - 1	27.8	9 - 1	25.4
11 - 1	26.8	11 - 1	27.5	11 - 1	25.7
12 - 1	27	12 - 1	27.9	12 - 1	25.5
19 - 1	26.7	19 - 1	27.5	19 - 1	26.1

Westgard	تفسیر	MCHC		آزمایش	NIHON KOHDEN سل کانترا	دستگاه	10	ماه	1402	سال
Normal Control	R&D	Mean	27.38	SD	0.33	Exp. Date	1402/12/29		Lot. No	
High Control	R&D	Mean	28.38	SD	0.51	Exp. Date	1402/12/29		Lot. No	
Low Control	R&D	Mean	25.61	SD	0.28	Exp. Date	1402/12/29		Lot. No	

