

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	2.25
SD	0.17
CV	7.37
Day-Run	Value
1 - 1	2.2
2 - 1	2.2
3 - 1	2.2
4 - 1	2
7 - 1	2.5
8 - 1	2.3
9 - 1	2.5
10 - 1	2
11 - 1	2.5
12 - 1	2.3
14 - 1	2.3
15 - 1	2.3
16 - 1	2.4
17 - 1	2.2
18 - 1	2.4
21 - 1	2.2
23 - 1	2.1
24 - 1	2.3
25 - 1	2
28 - 1	2
29 - 1	2.1
30 - 1	2.4

High - CV : 7.44	
Mean	3.43
SD	0.25
CV	7.44
Day-Run	Value
1 - 1	3.6
2 - 1	3.1
3 - 1	3.4
4 - 1	3.3
7 - 1	3.1
8 - 1	3.1
9 - 1	3.8
10 - 1	3
11 - 1	3.6
12 - 1	3.3
14 - 1	3.2
15 - 1	3.8
16 - 1	3.6
17 - 1	3.6
18 - 1	3.2
21 - 1	3.6
23 - 1	3.4
24 - 1	3.8
25 - 1	3.2
28 - 1	3.4
29 - 1	3.7
30 - 1	3.6

Westgard	تفسیر	Magnesium(Mg)	آزمایش	اتوآنالایزر BS-800	دستگاه	11	ماه	1402	سال
Normal Control	biorex	Mean	2.26	SD	0.13	Exp. Date	1403/1/17	Lot. No	035

