

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 6.04	
Mean	2.26	Mean	3.24
SD	0.13	SD	0.2
CV	5.57	CV	6.04
Day-Run	Value	Day-Run	Value
2 - 1	2.1	2 - 1	3
3 - 1	2.2	3 - 1	3
4 - 1	2.2	4 - 1	3.3
5 - 1	2.3	5 - 1	3.4
6 - 1	2.3	6 - 1	3.1
7 - 1	2.4	7 - 1	3
9 - 1	2.3	9 - 1	3.3
10 - 1	2.1	10 - 1	3.4
11 - 1	2.3	11 - 1	3.2
12 - 1	2.2	12 - 1	3.2
13 - 1	2.3	13 - 1	3
14 - 1	2.4	14 - 1	3
16 - 1	2.4	16 - 1	3.3
17 - 1	2.2	17 - 1	3
18 - 1	2.5	18 - 1	3.2
19 - 1	2.5	19 - 1	3.8
20 - 1	2.2	20 - 1	3.1
21 - 1	2	21 - 1	3.3
23 - 1	2.1	23 - 1	3.3
24 - 1	2.2	24 - 1	3.4
25 - 1	2.3	25 - 1	3.2
26 - 1	2.1	26 - 1	3.4
27 - 1	2.3	27 - 1	3.5
28 - 1	2.3	28 - 1	3.4
30 - 1	2.3	30 - 1	3.3

Westgard	تفسیر	Magnesium(Mg)	آزمایش	اتوآنالایزر BS-800	دستگاه	10	ماه	1402	سال
Normal Control	biorex	Mean	2.22	SD	0.16	Exp. Date	1403/1/17	Lot. No	035

