

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	148.1	
2 - 1	147.01	1:2S
3 - 1	148.3	
4 - 1	147.3	
7 - 1	147.9	
8 - 1	147.9	
9 - 1	146.8	
10 - 1	147.6	1:2S
12 - 1	148.3	
14 - 1	147.1	
15 - 1	147.5	
16 - 1	147.5	
17 - 1	148.5	
18 - 1	148	
21 - 1	147.6	1:2S
23 - 1	148.2	
24 - 1	148.2	
25 - 1	147.3	
28 - 1	147.8	
29 - 1	147.3	
30 - 1	147.5	1:2S

Westgard	تفسیر	Sodium (Na)	آزمایش	AUDICOM	دستگاه	11	ماه	1402	سال
Control	زیست شیمی	Mean	147.64	SD	0.56	Exp. Date	1404/1/14	Lot. No	2103061-4
Next Month	Mean	147.77	SD	0.54	CV	0.36			

