

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal		High		Low	
Mean	273.78	Mean	467.96	Mean	86.68
SD	9.59	SD	9.78	SD	3.67
CV	3.5	CV	2.09	CV	4.24
Day-Run	Value	Day-Run	Value	Day-Run	Value
23 - 1	277	23 - 1	473	23 - 1	86.4
24 - 1	263	24 - 1	459	24 - 1	88
25 - 1	268.7	25 - 1	451.2	25 - 1	83.2
27 - 1	265	26 - 1	467.5	27 - 1	88.5
28 - 1	284	27 - 1	470	28 - 1	82
30 - 1	285	28 - 1	477	30 - 1	92
		30 - 1	478		

Westgard	تفسیر	PLT	آزمایش	NIHON KOHDEN سل کانتر	دستگاه	10	ماه	1402	سال
Normal Control	R&D	Mean	272.7	SD	10.68	Exp. Date	1402/10/25	Lot. No	
High Control	R&D	Mean	469.7	SD	10.12	Exp. Date	1402/10/25	Lot. No	
Low Control	R&D	Mean	87.5	SD	4.7	Exp. Date	1402/10/25	Lot. No	

