

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	2.9	
3 - 1	3.1	
4 - 1	3	
5 - 1	3.15	
6 - 1	2.95	
7 - 1	2.99	
9 - 1	3.4	1:2S
10 - 1	2.9	
11 - 1	2.93	
12 - 1	3.12	
13 - 1	3.15	
14 - 1	2.98	
16 - 1	2.99	
17 - 1	3.5	1:2S
18 - 1	3.12	
19 - 1	3.15	
20 - 1	2.97	
21 - 1	2.99	
23 - 1	3.24	
24 - 1	3.26	
25 - 1	2.91	
26 - 1	3.1	
27 - 1	3.24	
28 - 1	2.98	
30 - 1	3.12	

Westgard	تفسیر	PSA	آزمایش	2000 ایمولایت XPI-	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	3.02	SD	0.18	Exp. Date	1403/8/30	Lot. No	365-02-028
Next Month	Mean	3.09	SD	0.16	CV	5.03			

