

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 2.36	
Mean	3.47	Mean	8.18
SD	0.1	SD	0.19
CV	2.84	CV	2.36
Day-Run	Value	Day-Run	Value
2 - 1	3.6	2 - 1	8.1
3 - 1	3.5	3 - 1	8.3
4 - 1	3.5	4 - 1	8.2
5 - 1	3.54	5 - 1	8
6 - 1	3.5	6 - 1	8.3
7 - 1	3.5	7 - 1	8
9 - 1	3.5	9 - 1	8.2
10 - 1	3.4	10 - 1	8.2
11 - 1	3.4	11 - 1	8.2
12 - 1	3.4	12 - 1	8
13 - 1	3.5	13 - 1	8.5
14 - 1	3.3	14 - 1	8
16 - 1	3.5	16 - 1	8
17 - 1	3.3	17 - 1	8.6
18 - 1	3.5	18 - 1	8.1
19 - 1	3.4	19 - 1	8.2
20 - 1	3.6	20 - 1	8.3
21 - 1	3.3	21 - 1	8
23 - 1	3.5	23 - 1	8.1
24 - 1	3.4	24 - 1	8.4
25 - 1	3.5	25 - 1	8.3
26 - 1	3.6	26 - 1	8.1
27 - 1	3.5	27 - 1	7.9
28 - 1	3.4	28 - 1	8.6
30 - 1	3.7	30 - 1	8

Westgard	تفسیر	Phosphorus (P)	آزمایش	اتوآنالایزر BS-800	دستگاه	10	ماه	1402	سال
Normal Control	biorex	Mean	3.5	SD	0.1	Exp. Date	1403/1/17	Lot. No	035

