

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	9.2	
2 - 1	10.62	
3 - 1	10.2	
4 - 1	10.7	
7 - 1	10.3	
8 - 1	9.75	
9 - 1	9.45	
10 - 1	10.51	
11 - 1	10.33	
12 - 1	9.2	
14 - 1	9.63	
15 - 1	9.25	
16 - 1	10.2	
17 - 1	9.85	
18 - 1	10.7	
21 - 1	11	1:2S
23 - 1	10.3	1:3S Cusum Reject
24 - 1	9.72	
25 - 1	9.2	
26 - 1	9.35	
28 - 1	9.98	
29 - 1	10.4	
30 - 1	10.2	

Westgard	تفسیر	Progesterone	آزمایش	SIEMENS Centaur CP	دستگاه	11	ماه	1402	سال
Control	biorex	Mean	9.86	SD	0.51	Exp. Date	1403/6/19	Lot. No	36-02-028
Next Month	Mean	10	SD	0.55	CV	5.52			

