

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	9.48	
3 - 1	10.32	
4 - 1	9.35	
5 - 1	9.22	
6 - 1	10.5	
7 - 1	10.12	
9 - 1	9.4	
10 - 1	9.65	
11 - 1	9.68	
12 - 1	10.2	
13 - 1	9.8	
14 - 1	10.5	
16 - 1	9.57	
17 - 1	9.2	
18 - 1	9.45	
19 - 1	10.7	
20 - 1	10.6	
21 - 1	10.7	
23 - 1	9.8	4:1S
24 - 1	9.57	
25 - 1	10.23	
26 - 1	10.1	
27 - 1	9.35	
28 - 1	9.2	
30 - 1	9.42	

Westgard	تفسیر	Progesterone	آزمایش	SIEMENS Centaur CP	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	9.44	SD	0.7	Exp. Date	1403/6/19	Lot. No	36-02-028
Next Month	Mean	9.86	SD	0.51	CV	5.19			

