

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	19.4	
2 - 1	19.1	
3 - 1	18.4	
4 - 1	19.9	1:2S
7 - 1	17.9	
8 - 1	17.9	
9 - 1	19.3	
10 - 1	19.4	
11 - 1	17.9	
12 - 1	18.5	
14 - 1	19.1	
15 - 1	19.3	
16 - 1	18.8	
17 - 1	18.2	
18 - 1	18.1	
21 - 1	18.9	
23 - 1	19.7	1:2S
24 - 1	19.1	
25 - 1	19.2	
26 - 2	17.9	
28 - 1	18.9	
29 - 1	18.1	
30 - 1	19.1	

Westgard	تفسیر	Prolactin	آزمایش	2000 ایمولایت XPI-	دستگاه	11	ماه	1402	سال
Control	biorex	Mean	18.45	SD	0.53	Exp. Date	1403/2/10	Lot. No	365-02-028
Next Month	Mean	18.79	SD	0.63	CV	3.33			

