

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	17.8	
3 - 1	18.3	
4 - 1	18.5	
5 - 1	17.9	
6 - 1	17.4	
7 - 1	18.9	
9 - 1	19.1	
10 - 1	18.1	
11 - 1	18.6	
12 - 1	19.2	
13 - 1	17.9	
14 - 1	17.8	
16 - 1	18.6	
17 - 1	18.8	
18 - 1	19.4	
19 - 1	19.2	
20 - 1	18.5	
21 - 1	18.2	
23 - 1	18.1	
24 - 1	18.5	
25 - 1	19.1	
26 - 1	18.4	
27 - 1	18.2	
28 - 1	17.9	
30 - 1	18.8	

Westgard	تفسیر	Prolactin	آزمایش	2000 ایمولایت XPI-	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	18.17	SD	0.99	Exp. Date	1403/2/10	Lot. No	365-02-028
Next Month	Mean	18.45	SD	0.53	CV	2.85			

