

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	102.8	
3 - 1	103.9	
4 - 1	107.3	
5 - 1	108	
6 - 1	106.9	
7 - 1	107.3	
9 - 1	105.9	
10 - 1	108	
11 - 1	109.1	
12 - 1	106.8	
13 - 1	103.2	
14 - 1	109.2	
16 - 1	103.7	
17 - 1	105.3	
18 - 1	109	
19 - 1	105.9	1:3S
20 - 1	105.2	
21 - 1	109.4	
23 - 1	108.1	
24 - 1	106.2	
25 - 1	107.3	
26 - 1	106.8	
27 - 1	109.5	
28 - 1	110.3	
30 - 1	107.02	

Westgard	تفسیر	RF	آزمایش	BS-480 اتوآنالایزر	دستگاه	10	ماه	1402	سال
Control	aptec	Mean	107.07	SD	2.31	Exp. Date	1403/8/19	Lot. No	06095
Next Month	Mean	106.88	SD	2.05	CV	1.92			

