

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

| Normal  |       |
|---------|-------|
| Mean    | 1.01  |
| SD      | 0.08  |
| CV      | 8.02  |
| Day-Run | Value |
| 1 - 1   | 0.92  |
| 2 - 1   | 1.15  |
| 3 - 1   | 1.05  |
| 4 - 1   | 1.04  |
| 7 - 1   | 1.05  |
| 8 - 1   | 1     |
| 9 - 1   | 0.95  |
| 10 - 1  | 0.92  |
| 11 - 1  | 1.05  |
| 12 - 1  | 1.1   |
| 14 - 1  | 0.93  |
| 15 - 1  | 1.1   |
| 15 - 1  | 1.1   |
| 16 - 1  | 0.9   |
| 17 - 1  | 1.05  |
| 18 - 1  | 0.95  |
| 18 - 1  | 0.95  |
| 21 - 1  | 1.07  |
| 23 - 1  | 1.1   |
| 24 - 1  | 0.9   |
| 25 - 1  | 1.1   |
| 28 - 1  | 0.95  |
| 29 - 1  | 1.1   |
| 30 - 1  | 0.9   |

| High - CV : 2.56 |       |
|------------------|-------|
| Mean             | 4.26  |
| SD               | 0.11  |
| CV               | 2.56  |
| Day-Run          | Value |
| 1 - 1            | 4.44  |
| 2 - 1            | 4.3   |
| 3 - 1            | 4.11  |
| 4 - 1            | 4.32  |
| 7 - 1            | 4.2   |
| 8 - 1            | 4.25  |
| 9 - 1            | 4.2   |
| 10 - 1           | 4.18  |
| 11 - 1           | 4.4   |
| 12 - 1           | 4.15  |
| 14 - 1           | 4.36  |
| 15 - 1           | 4.3   |
| 16 - 1           | 4.1   |
| 17 - 1           | 4.2   |
| 17 - 2           | 4.2   |
| 18 - 1           | 4.2   |
| 21 - 1           | 4.5   |
| 23 - 1           | 4.33  |
| 24 - 1           | 4.15  |
| 25 - 1           | 4.15  |
| 28 - 1           | 4.26  |
| 29 - 1           | 4.35  |
| 30 - 1           | 4.38  |

|                |        |        |        |                    |        |           |           |         |     |
|----------------|--------|--------|--------|--------------------|--------|-----------|-----------|---------|-----|
| Westgard       | تفسیر  | Bili T | آزمایش | اتوآنالایزر BS-800 | دستگاه | 11        | ماه       | 1402    | سال |
| Normal Control | biorex | Mean   | 1.02   | SD                 | 0.06   | Exp. Date | 1403/1/17 | Lot. No | 035 |

