

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	7.08
SD	0.1
CV	1.35
Day-Run	Value
1 - 1	7
2 - 1	7.1
3 - 1	7
4 - 1	7
7 - 1	7
8 - 1	7
9 - 1	7
10 - 1	7
11 - 1	7
12 - 1	7
14 - 1	7.1
15 - 1	7.1
16 - 1	7.1
17 - 1	7.2
18 - 1	7.2
21 - 1	7.2
23 - 1	6.9
24 - 1	7.2
25 - 1	7.1
28 - 1	7.2
29 - 1	7.2
30 - 1	7.2

High - CV : 2.98	
Mean	8.17
SD	0.24
CV	2.98
Day-Run	Value
1 - 1	8.5
2 - 1	8.4
3 - 1	8.3
4 - 1	8.4
7 - 1	8.2
8 - 1	8.4
9 - 1	7.9
10 - 1	8.1
11 - 1	8.4
12 - 1	8
14 - 1	8
15 - 1	8
16 - 1	8.1
17 - 1	7.9
18 - 1	7.8
21 - 1	8
23 - 1	8.1
24 - 1	7.9
25 - 1	8
28 - 1	8.2
29 - 1	8.5
30 - 1	8.7

Westgard	تفسیر	Protein Total	آزمایش	اتوآنالایزر BS-800	دستگاه	11	ماه	1402	سال
Normal Control	biorex	Mean	7	SD	0.2	Exp. Date	1403/1/17	Lot. No	035

