

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	2.43	
2 - 1	2.36	
3 - 1	2.1	
4 - 1	2.22	
7 - 1	2.1	
8 - 1	2.3	
9 - 1	2.45	
10 - 1	2.5	
11 - 1	2.75	
12 - 1	2.6	
14 - 1	2.2	
15 - 1	2.25	
16 - 1	2.5	
17 - 1	2.1	
18 - 1	2	1:2S
21 - 1	2.6	
23 - 1	2.7	
24 - 1	2.25	
25 - 1	2.24	
26 - 1	2.36	
28 - 1	2.3	
29 - 1	2.45	
30 - 1	2.65	

Westgard	تفسیر	T3	آزمایش	SIEMENS Centaur CP	دستگاه	11	ماه	1402	سال
Control	biorex	Mean	2.46	SD	0.21	Exp. Date	1403/6/19	Lot. No	36-02-028
Next Month	Mean	2.37	SD	0.21	CV	8.78			

