

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	2.99	1:2S
3 - 1	2.34	
4 - 1	2.27	
5 - 1	2.54	
6 - 1	2.76	
7 - 1	2.32	
9 - 1	2.34	
10 - 1	2.56	
11 - 1	2.42	
12 - 1	2.45	
13 - 1	2.75	
14 - 1	2.8	
16 - 1	2.25	
17 - 1	2.3	
18 - 1	2.1	
19 - 1	2.35	
20 - 1	2.24	
21 - 1	2.44	
23 - 1	2.56	
24 - 1	2.43	
25 - 1	2.77	
26 - 1	2.5	
27 - 1	2.48	
28 - 1	2.35	
30 - 1	2.27	

Westgard	تفسیر	T3	آزمایش	SIEMENS Centaur CP	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	2.43	SD	0.19	Exp. Date	1403/6/19	Lot. No	36-02-028
Next Month	Mean	2.46	SD	0.21	CV	8.65			

