

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal		High - CV : 1.16	
Mean	172.01	Mean	231.36
SD	2.57	SD	2.69
CV	1.49	CV	1.16
Day-Run	Value	Day-Run	Value
2 - 1	170	2 - 1	230
3 - 1	173	3 - 1	235
4 - 1	172	4 - 1	230
5 - 1	170	5 - 1	233
6 - 1	172	6 - 1	233
7 - 1	173	7 - 1	230
9 - 1	169	9 - 1	231
10 - 1	170	10 - 1	233
11 - 1	172	11 - 1	236
12 - 1	171	12 - 1	232
13 - 1	170	13 - 1	234
14 - 1	170	14 - 1	233
17 - 1	174	17 - 1	228
17 - 2	175	18 - 1	235
18 - 1	171	19 - 1	229
19 - 1	172	20 - 1	232
20 - 1	177	21 - 1	230
23 - 1	178	24 - 1	226
24 - 1	170	25 - 1	230.2
26 - 1	174.2	26 - 1	232
27 - 1	175	27 - 1	234
28 - 1	170	28 - 1	227
30 - 1	168	30 - 1	228

Westgard	تفسیر	Triglycerides	آزمایش	اتوآنالایزر BS-800	دستگاه	10	ماه	1402	سال
Normal Control	biorex	Mean	173.21	SD	2.64	Exp. Date	1403/1/17	Lot. No	035

