

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 1	311	
3 - 1	320	1:3S
4 - 1	319	1:3S
5 - 1	310.5	
6 - 1	305.1	
7 - 1	310	
9 - 1	317	
10 - 1	308.5	
11 - 1	307.5	
12 - 1	317.5	
13 - 1	312	1:3S
14 - 1	304	
16 - 1	316	
17 - 1	308	
18 - 1	300.5	1:2S
19 - 1	303.6	
20 - 1	301.8	
21 - 1	307	
24 - 1	315	
25 - 1	301	
26 - 1	311	
27 - 1	312	
28 - 1	310.6	
30 - 1	315.4	

Westgard	تفسیر	TIBC	آزمایش	اتوآنالایزر BS-800	دستگاه	10	ماه	1402	سال
Control	biorex	Mean	310	SD	4.5	Exp. Date	1403/1/17	Lot. No	035
Next Month	Mean	310.17	SD	5.7	CV	1.84			

