

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
1 - 1	521	
2 - 1	526	
3 - 1	543	
4 - 1	512	
7 - 1	500	
8 - 1	512	
9 - 1	495	
10 - 1	490	
11 - 1	521	
12 - 1	550	
14 - 1	525	
15 - 1	520	
16 - 1	524	
17 - 1	501	
18 - 1	508	
21 - 1	495	
23 - 1	490	
24 - 1	506	
25 - 1	510	
26 - 2	523	
28 - 1	502	
29 - 1	513	
30 - 1	521	

Westgard	تفسیر	Testosterone	آزمایش	2000 ایمولایت XPI-	دستگاه	11	ماه	1402	سال
Control	Labex	Mean	512.36	SD	22.44	Exp. Date	1402/11/15	Lot. No	36502028
Next Month	Mean	513.39	SD	15.47	CV	3.01			

