

آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

Control Serum - Levey Jennings Chart

Normal	
Mean	34.17
SD	0.96
CV	2.81
Day-Run	Value
1 - 1	35.4
2 - 1	35
3 - 1	33.1
4 - 1	34
7 - 1	33.1
8 - 1	33
9 - 1	33.5
10 - 1	35
11 - 1	34.6
12 - 1	33
14 - 1	35.2
15 - 1	33
16 - 1	34.4
17 - 1	33.5
18 - 1	35
21 - 1	35.1
23 - 1	35.5
24 - 1	32.8
25 - 1	33.3
28 - 1	34.1
29 - 1	35.4
30 - 1	34.8

High - CV : 1.06	
Mean	130.81
SD	1.39
CV	1.06
Day-Run	Value
1 - 1	133.2
2 - 1	130
3 - 1	128.1
4 - 1	131
7 - 1	132
8 - 1	128
9 - 1	130
10 - 1	130
11 - 1	132
12 - 1	131.2
14 - 1	131
15 - 1	132
16 - 1	132
17 - 1	130
18 - 1	130
21 - 1	132
23 - 1	130
24 - 1	132
25 - 1	129
28 - 1	132
29 - 1	130.1
30 - 1	132.2

Westgard	تفسیر	Urea	آزمایش	اتوآنالایزر BS-800	دستگاه	11	ماه	1402	سال
Normal Control	biorex	Mean	33.79	SD	1	Exp. Date	1403/1/17	Lot. No	035

