

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Levey Jennings Chart

Normal	
Mean	33.79
SD	0.85
CV	2.5
Day-Run	Value
2 - 1	34.1
3 - 1	34.6
4 - 1	34
5 - 1	34.3
6 - 1	33
7 - 1	35
9 - 1	34
10 - 1	34.4
11 - 1	33.6
12 - 1	34
13 - 1	34.5
14 - 1	32.3
16 - 1	33
17 - 1	35
18 - 1	34.5
19 - 1	33
20 - 1	33
21 - 1	34.7
24 - 1	35
25 - 1	33.3
26 - 1	33.2
27 - 1	33.2
28 - 1	33.1
30 - 1	32.2

High - CV : 1.39	
Mean	131.33
SD	1.83
CV	1.39
Day-Run	Value
2 - 1	130
3 - 1	134
4 - 1	133
5 - 1	133.1
6 - 1	130
7 - 1	128
9 - 1	130
10 - 1	131
11 - 1	133.7
12 - 1	130.6
13 - 1	130.1
14 - 1	129
16 - 1	131.9
17 - 1	130.9
18 - 1	135
19 - 1	130
20 - 1	134.1
21 - 1	132.6
24 - 1	131
25 - 1	130
26 - 1	132
27 - 1	128.8
28 - 1	131.9
30 - 1	131.3

Westgard	تفسیر	Urea	آزمایش	اتوآنالایزر BS-800	دستگاه	10	ماه	1402	سال
Normal Control	biorex	Mean	33.44	SD	0.9	Exp. Date	1403/1/17	Lot. No	035

