

# آزمایشگاه پاتوبیولوژی و ژنتیک کاوش گرگان

## Control Serum - Normal Levey Jennings Chart

Day-Run	Value	Description
2 - 2	550	
3 - 1	556	
4 - 1	578	
7 - 1	574	
8 - 1	570	
9 - 1	581	
10 - 1	588	
11 - 1	582	
12 - 1	590	
14 - 1	582	
15 - 1	578	
17 - 1	560	
18 - 1	552	
21 - 1	565	
23 - 1	593	
24 - 1	575	
25 - 1	586	
26 - 1	590	
28 - 1	586	
30 - 1	570	

Westgard	تفسیر	Vitamin B12( Cobalamin )	آزمایش	Maglumi	دستگاه	11	ماه	1402	سال
Control		Maglumi	Mean	572.09	SD	13.21	Exp. Date	1402/11/25	Lot. No
Next Month		Mean	575.3	SD	13.03	CV	2.27		

